


LOCKING CANNULATED BLADE TRAY LAYOUT




T15/T20 Screw Instrumentation Base Tray

Key	Description	Qty	Item#
1	2.5mm dia. Neutral & Load End Green/Gold Drill Guide	1	01-1200-0054
2	3.5mm/2.5mm Double Drill Guide	1	01-1050-0009
3	2.5mm Threaded Drill Guide	3	01-1200-0067
4	Small Bone Clamp	1	01-1200-0057
4	Infant Bone Clamp	1	01-1200-0074
5	2.0 X 150mm Guide Wire	4	01-0907-0022 Cobalt Chrome
-	-	-	01-9473-5115 Stainless Steel
6	3.5mm Depth Gauge Sleeve (long)	1	01-1200-0078
7	Mini In-Line Ratchet	2	01-1030-001
8	Screw Forceps	1	01-1030-007
9	Precision Wire Guide- 2.0mm	1	01-1200-0014
10	2.5 Calibrated Drill Bit	2	01-1050-0032
11	T15 Hexalobe Retaining Driver, Short	1	01-0903-0005
11	T15 Hexalobe Retaining Driver, Long	1	01-1200-0087
12	Drill Bit 2.5mm	2	01-1050-0002
14	Triangular Positioning Plate 90-40-50 Degrees	1	01-1200-0069
15	Triangular Positioning Plate 80-70-30 Degrees	1	01-1200-0070
16	Triangular Positioning Plate 100-60-20 Degrees	1	01-1200-0071

NOTE: The stainless steel guide wire options listed may not be available in all countries.

Key	Description	Qty	Item#
17	4.5mm dia. Neutral & Load End Green/Gold Drill Guide	1	01-1200-0055
18	Double Drill Sleeve 4.5mm/3.2mm	1	01-1200-0056
19	3.2mm Threaded Drill Guide	3	01-1200-0042
20	3.2mm Drill Bit, Calibrated	2	01-1200-0041
20	3.2mm Drill Bit	2	01-1200-0051
21	T20 Hexalobe Retaining Driver, Short	1	01-1200-0088
21	T20 Hexalobe Retaining Driver, Long	1	01-1200-0089
22	Cortical Tap 3.5mm	1	01-1050-0006
23	4.5mm Tap	1	01-1200-0052
24	Large Bone Clamp	1	01-1200-0058
25	4.5mm Depth Gauge Sleeve (long)	1	01-1200-0079
26	Long insert (0-100mm)	1	01-0999-2001
27	Cleaning Stylet	1	01-1010-004
28	Bending Iron-Right	1	01-1200-0062
29	Bending Iron-Left	1	01-1200-0064


LOCKING CANNULATED BLADE TRAY LAYOUT



Locking Cannulated Blade Plate (LCB) Instrumentation Base Tray

Key	Description	Qty	Item#
1	3.5mm Cannulated Chisel-Calibrated	1	01-1200-0007
2	2.7mm Cannulated Chisel-Calibrated	1	01-1200-0006
3	Spatula Spacer, Child/Infant	1	01-1200-0075
4	2.7mm Threaded Inserter	1	01-1200-0009
5	Chisel Guide, Child	1	01-1200-0002
6	Precision Wire Guide- 2.4mm	1	01-1200-0115
7	2.4mm Threaded Guide Wire 150mm	4	01-1200-0050
8	4.5mm Cannulated Chisel-Calibrated	1	01-1200-0008
9	Spatula Spacer, Adolescent	1	01-1200-0060
10	4.5mm Threaded Inserter	1	01-1200-0011
11	Slap Hammer, Small	1	01-1200-0076
12	Chisel Guide, Adolescent	1	01-1200-0003
13	Tuning Fork	1	01-1200-0026


LOCKING CANNULATED BLADE TRAY LAYOUT



90° and 100° Locking Cannulated Blade Plates (LCB), Infant and Child (3.5mm)

Key	Description (Blade Length X Offset X Number of Holes)	Qty	Item#
1	90 Deg. Cann Infant Blade Plate 25mm x 5mm x 3 Hole	2	00-1200-1000
2	90 Deg. Cann Infant Blade Plate 30mm x 5mm x 3 Hole	2	00-1200-1001
3	90 Deg. Cann Infant Blade Plate 35mm x 5mm x 3 Hole	2	00-1200-1002
4	90° Locking Cannulated Child Blade 25mm x 6mm x 3 Hole	2	00-1200-3500
5	90° Locking Cannulated Child Blade 30mm x 6mm x 3 Hole	2	00-1200-3501
6	90° Locking Cannulated Child Blade 35mm x 6mm x 3 Hole	2	00-1200-3502
7	90 °Locking Cannulated Child Blade 40mm x 6mm x 3 Hole	2	00-1200-3503
8	90° Locking Cannulated Child Blade 45mm x 6mm x 3 Hole	2	00-1200-3504
9	90° Locking Cannulated Child Blade 50mm x 6mm x 3 Hole	2	00-1200-3505
10	90° Locking Cannulated Child Blade 25mm x 10mm x 3 Hole	2	00-1200-3506
11	90° Locking Cannulated Child Blade 30mm x 10mm x 3 Hole	2	00-1200-3507
12	90° Locking Cannulated Child Blade 35mm x 10mm x 3 Hole	2	00-1200-3508
13	90° Locking Cannulated Child Blade 40mm x 10mm x 3 Hole	2	00-1200-3509
14	90° Locking Cannulated Child Blade 45mm x 10mm x 3 Hole	2	00-1200-3510
15	90° Locking Cannulated Child Blade 50mm x 10mm x 3 Hole	2	00-1200-3511
16	100° Locking Cannulated Child Blade 25mm x 6mm x 3 Hole	2	00-1200-4000
17	100° Locking Cannulated Child Blade 30mm x 6mm x 3 Hole	2	00-1200-4001


Key	Description	Qty	Item#
18	100° Locking Cannulated Child Blade 35mm x 6mm x 3 Hole	2	00-1200-4002
19	100° Locking Cannulated Child Blade 40mm x 6mm x 3 Hole	2	00-1200-4003
20	100° Locking Cannulated Child Blade 45mm x 6mm x 3 Hole	2	00-1200-4004
21	100° Locking Cannulated Child Blade 50mm x 6mm x 3 Hole	2	00-1200-4005



90° Locking Cannulated Blade Plates (LCB), Adolescent (4.5mm)

Key	Description (Blade Length X Offset X Number of Holes)	Qty	Item#
1	90° Locking Cannulated Adolescent Blade 40mm x 6mm x 4 Hole	2	00-1200-4500
2	90° Locking Cannulated Adolescent Blade 45mm x 6mm x 4 Hole	2	00-1200-4501
3	90° Locking Cannulated Adolescent Blade 50mm x 6mm x 4 Hole	2	00-1200-4502
4	90° Locking Cannulated Adolescent Blade 55mm x 6mm x 4 Hole	2	00-1200-4503
5	90° Locking Cannulated Adolescent Blade 60mm x 6mm x 4 Hole	2	00-1200-4504
6	90° Locking Cannulated Adolescent Blade 65mm x 6mm x 4 Hole	2	00-1200-4505
7	90° Locking Cannulated Adolescent Blade 70mm x 6mm x 4 Hole	2	00-1200-4506
8	90° Locking Cannulated Adolescent Blade 40mm x 14mm x 4 Hole	2	00-1200-4514
9	90° Locking Cannulated Adolescent Blade 45mm x 14mm x 4 Hole	2	00-1200-4515
10	90° Locking Cannulated Adolescent Blade 50mm x 14mm x 4 Hole	2	00-1200-4516
11	90° Locking Cannulated Adolescent Blade 55mm x 14mm x 4 Hole	2	00-1200-4517
12	90° Locking Cannulated Adolescent Blade 60mm x 14mm x 4 Hole	2	00-1200-4518
13	90° Locking Cannulated Adolescent Blade 65mm x 14mm x 4 Hole	2	00-1200-4519
14	90° Locking Cannulated Adolescent Blade 70mm x 14mm x 4 Hole	2	00-1200-4520


LOCKING CANNULATED BLADE TRAY LAYOUT



130° Locking Cannulated Blade Plates (LCB), Child and Adolescent (3.5mm and 4.5mm)

Key	Description (Blade Length X Offset X Number of Holes)	Qty	Item#
1	130° Locking Cannulated Child Blade 40mm x 0mm x 4 Hole	2	00-1200-5000
2	130° Locking Cannulated Child Blade 45mm x 0mm x 4 Hole	2	00-1200-5001
3	130° Locking Cannulated Child Blade 50mm x 0mm x 4 Hole	2	00-1200-5002
4	130° Locking Cannulated Child Blade 55mm x 0mm x 4 Hole	2	00-1200-5003
5	130° Locking Cannulated Child Blade 60mm x 0mm x 4 Hole	2	00-1200-5004
6	130° Locking Cannulated Adolescent Blade 45mm x 0mm x 4 Hole	2	00-1200-7000
7	130° Locking Cannulated Adolescent Blade 50mm x 0mm x 4 Hole	2	00-1200-7001
8	130° Locking Cannulated Adolescent Blade 55mm x 0mm x 4 Hole	2	00-1200-7002
9	130° Locking Cannulated Adolescent Blade 60mm x 0mm x 4 Hole	2	00-1200-7003
10	130° Locking Cannulated Adolescent Blade 65mm x 0mm x 4 Hole	2	00-1200-7004
11	130° Locking Cannulated Adolescent Blade 70mm x 0mm x 4 Hole	2	00-1200-7005
12	130° Locking Cannulated Adolescent Blade 75mm x 0mm x 4 Hole	2	00-1200-7006
13	130° Locking Cannulated Adolescent Blade 80mm x 0mm x 4 Hole	2	00-1200-7007

LOCKING CANNULATED BLADE TRAY LAYOUT




3.5mm Bilateral Modular Screw Caddy

Key	Description	Qty	Item#
1	3.5mm Non-Locking Cortical Screw	-	00-0903-2510
	10mm - 30mm in 2mm increments	6	00-0903-2530
	3.5mm Non-Locking Cortical Screw	-	00-0903-2532
	32mm - 50mm in 2mm increments	4	00-0903-2550
	3.5mm Non-Locking Cortical Screw	-	00-0903-2555
	55mm - 70mm in 5mm increments	4	00-0903-2570

Key	Description	Qty	Item#
2	3.5mm Locking Cortical Screw	-	00-0903-2610
	10mm - 30mm in 2mm increments	6	00-0903-2630
	3.5mm Locking Cortical Screw	-	00-0903-2632
	32mm - 50mm in 2mm increments	4	00-0903-2650
	3.5mm Locking Cortical Screw	-	00-0903-2655
	55mm - 70mm in 5mm increments	4	00-0903-2670

Key	Description	Qty	Item#
3	3.5mm Cannulated Locking Cortical Screw	-	00-0907-3730
	30mm - 40mm in 2mm increments	4	00-0907-3740
	3.5mm Cannulated Locking Cortical Screw	-	00-0907-3745
	45mm - 70mm in 5mm increments	4	00-0907-3770

LOCKING CANNULATED BLADE TRAY LAYOUT



4.5mm Bilateral Modular Screw Caddy

Key	Description	Qty	Item#
1	4.5mm Non-Locking Cortical Screw	-	00-0907-4510
	10mm - 50mm in 2mm increments	5	00-0907-4550
	4.5mm Non-Locking Cortical Screw	-	00-0907-4555
	50mm - 80mm in 5mm increments	4	00-0907-4580

Key	Description	Qty	Item#
2	4.5mm Locking Cortical Screw	-	00-0907-4610
	10mm - 40mm in 2mm increments	5	00-0907-4640
	4.5mm Locking Cortical Screw	-	00-0907-4642
	42mm - 50mm in 2mm increments	4	00-0907-4650
	4.5mm Locking Cortical Screw	-	00-0907-4655
	55mm - 80mm in 5mm increments	4	00-0907-4680

Key	Description	Qty	Item#
3	4.5mm Cannulated Locking Cortical Screw	-	00-0907-4730
	30mm - 40mm in 2mm increments	4	00-0907-4740
	4.5mm Cannulated Locking Cortical Screw	-	00-0907-4745
	45mm - 80mm in 5mm increments	4	00-0907-4780

